


# VETERINARY SCREWS, DRILL BITS AND TAPS.

# Standard Screws, Drill Bits and Taps for Veterinary Use



SCREW DIAMETER (mm)												
Thread Diameter	1.5	2.0	2.4	2.7	3.5	4.5	5.5	4.0			6.5	
Screw Type	Cortex							Cancellous				
								Full thread	Partial thread	Full thread	16 mm thread	32 mm thread
Vet Part Number	VS102	VS202	VS203	VS205	VS302	VS402	–	VS403	VS404	VS601	VS602	VS603
Self-tapping	VS101	VS201	–	VS204	VS301	VS401	VS501					
Non-self-tapping												
DRILL BIT and TAP DIAMETER (mm)												
Drill Bit for Gliding Hole	1.5	2.0	2.4	2.7	3.5	4.5	5.5	–			4.5	
	⊗	⊗	⊗	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙
Drill Bit for Threaded Hole	1.1	1.5	1.8	2.0	2.5	3.2	4.0	2.5			3.2	
Tap Catalog Number	1.5 311.15 –	2.0 311.19 311.21 (long)	–	2.7 311.26 311.28 (long)	3.5 311.32 311.33 (long)	4.5 311.46 311.48	5.5 311.55	4.0 311.34			6.5 311.66	
Drive Type	1.5 mm/2.0 mm Cruciform		2.4 mm Cruciform	2.5 mm Hex		3.5 mm Hex		2.5 mm Hex			3.5 mm Hex	

# Locking Screws and Drill Bits for Veterinary Use



SCREW DIAMETER (mm)							
Thread Diameter	1.5	2.0	2.4	2.7	3.5	4.0	5.0
Screw Type	Locking Head						
Vet Part Number	N/A	VS207	VS208	VS206	VS303	VS406	VS502
Self-tapping							
DRILL BIT DIAMETER (mm)							
	⊙	⊙	⊙	⊙	⊙	⊙	⊙
Drill Bit for Threaded Hole	1.1	1.5	1.8	2.0	2.8	3.2	4.3
Drive Type	T4 StarDrive	T6 StarDrive	T8 StarDrive		T15 StarDrive	T25 StarDrive	

Illustrations not actual size



COMPANIES OF Johnson & Johnson

**DePuy Synthes Vet**  
1302 Wrights Lane East  
West Chester, PA 19380  
Telephone: (610) 719-5000  
To order: (800) 523-0322

[www.depuyssynthes.com](http://www.depuyssynthes.com)

**Synthes (Canada) Ltd.**  
2566 Meadowpine Boulevard  
Mississauga, Ontario L5N 6P9  
Telephone: (905) 567-0440  
To order: (800) 668-1119  
Fax: (905) 567-3185